Effect of 7.0 Tesla MRI on Upper Eyelid Implants

Schrom, Thomas M.D.^{*}; Thelen, Ariane M.D.^{*}; Asbach, Patrick M.D.[†]; Bauknecht, Hans-Christian M.D.[†]

Author Information

*Department of Otorhinolaryngology and †Department of Radiology, Charité–Universitaetsmedizin Berlin, Charité Campus Mitte, Berlin, Germany

Accepted for publication March 31, 2006.

Address correspondence and reprint requests to Thomas Schrom, M.D., Department of Otorhinolaryngology, Charité -University Medicine Berlin, Campus Mitte, Schumannstrasse 20/21, 10117 Berlin, Germany; E-mail: thomas.schrom@gmx.de.

Ophthalmic Plastic & Reconstructive Surgery: November-December 2006 - Volume 22 - Issue 6 - p 480-482

doi: 10.1097/01.iop.0000240807.42320.d9

• Buy

Metrics	
Abstract	In Brief

To ensure the MRI compatibility of various eyelid implants in high-field MRI, 3 eyelid weights made of pure gold (99.99%), pure platinum (99.95%), and a platinum (97%)/iridium (3%) alloy were examined in vitro. Temperature changes, position changes, and imaging artifacts of the different implants were determined in a small-bore 7.0 Tesla MRI system. The 7.0 Tesla MRI system demonstrated that none of the eyelid implants carried a risk of heating or dislocation; therefore, these implants are MRI compatible up to a magnetic field strength of 7.0 Tesla.

Eyelid weights made of pure gold (99.99%), pure platinum (99.95%), or a platinum (97%)/iridium (3%) alloy are MR compatible up to a magnetic field strength of 7.0 Tesla.

©2006The American Society of Opthalmic Plastic and Reconstructive Surgery, Inc.

∧Back to Top



Never Miss an Issue

Get new journal Tables of Contents sent right to your email inbox Get New Issue Alerts